

# PHERASTAR<sup>®</sup> FSX

The gold standard for high-throughput screening



  
**BMG LABTECH**  
*The Microplate Reader Company*

[www.bmglabtech.com](http://www.bmglabtech.com)



## Your benefits at a glance:

- Best sensitivity and speed in all assays
- Assay Optic Modules for easiest optical configuration
- High temporal resolution in TR-FRET detection
- Cross-talk reduction for AlphaScreen and luminescence
- Steady temperature unaffected by external changes
- Automation friendly
- Optimisation tool for TRF/TR-FRET and AlphaScreen
- Made-in-Germany dependability

## The best performance in all assays

Developed to fulfil the hardest High-Throughput Screening (HTS) requirements, the PHERAstar® FSX high-end microplate reader meets the needs of every modern laboratory. Specifically conceived for the highest sensitivity and the fastest speed, its unique features make it superior to any other microplate reader currently on the market.

High-end German engineering combined with the latest technologies makes the PHERAstar FSX the gold standard microplate reader for HTS.

This multi-mode plate reader performs all the leading detection technologies:

- UV/vis absorbance
- Fluorescence intensity, including FRET
- Luminescence (flash and glow), including BRET
- Fluorescence polarization/anisotropy
- Time-resolved fluorescence (TRF), including TR-FRET
- AlphaScreen®, AlphaLISA® and AlphaPlex™

Whatever your application, the PHERAstar FSX provides uncompromised sensitivity, speed and dynamic range in all plate formats up to 3456 wells. Three dedicated light sources, Simultaneous Dual Emission, Decay Curve Monitoring and assay-dedicated Optic Modules are just some of its key features. Assay flexibility is enhanced by software-controlled top/bottom reading, on-board reagent injectors, precise temperature control and multi-mode shaking capabilities.

## Sensitivity and speed

The PHERAstar FSX combines the fast read times necessary for HTS with the sensitivity to read small volumes. The user can always find the best combination of sensitivity and speed by choosing the number of excitation flashes. In single flash mode, the reader can measure a 1536-well plate in 27 seconds - one of the fastest readers on the market. Even at low concentrations and small assay volumes, the unsurpassed sensitivity of its detection system provides outstanding signal-to-noise, %CV, and Z' values.

## Innovative technology

The outstanding sensitivity of the PHERAstar FSX is centered on an innovative optical design composed of a free air optical path, three independent light sources, Simultaneous Dual Emission detection, and high transmission filters. The optical system is able to adjust the detection focal height for both top and bottom readings to ensure the highest possible sample signal. Top or bottom detection can be directly selected by a mouse click, requiring no hardware changes. Depending on the application, users can choose one of the following light sources:

- High energy xenon flash lamp
- Laser for TRF/TR-FRET
- Laser for AlphaScreen, AlphaLISA and AlphaPlex

Four photomultiplier tubes, combined in two application dedicated matched pairs, are used as detectors. One pair is dedicated to simultaneous dual fluorescence and luminescence detection, the second to TRF-based signals.

## Switch between assays in seconds

Never worry about which filters or dichroic mirrors are installed. Our application-specific Optic Modules make the PHERAstar FSX the easiest reader to optically configure on the market. Optic Modules contain all the components required for a specific assay such as excitation and emission filters, dichroic mirrors, beam splitters and polarization filters. The modules are easy to use: they are automatically recognized via barcode and selected by the reader for the appropriate assay. The reader accommodates up to six Optic Modules. You can easily add or replace them within seconds.

## Cut read times in half

Numerous assays such as FRET, TR-FRET, BRET, FP and AlphaPlex™ require detection of two emission wavelengths. BMG LABTECH pioneered the technique of Simultaneous Dual Emission (SDE) detection for microplate readers. Thanks to SDE detection, the PHERAstar FSX can simultaneously detect two separate emission wavelengths in one single measurement. SDE detection can be used in any assay that measures two emission wavelengths or polarization vectors. It offers a significant speed advantage

by cutting read times in half and eliminates the typical drawbacks of double sequential detection such as flash-to-flash variations, photobleaching and decaying kinetic signals. SDE also reduces the variability caused by fluctuations in pH or temperature, and by evaporation.

## Ultra-fast detection of absorbance spectra


Why would you detect only a single wavelength if you could acquire a full spectrum in the same amount of time? Full-spectrum detection improves many absorbance assays as it can highlight shifting peaks or the presence of contaminants. Our spectrometer captures a full absorbance spectrum (220 - 1000 nm) at resolutions selectable from 1 to 10 nm in less than 1 second/well. Alternatively, up to 8 discrete wavelengths can be acquired instantaneously and simultaneously. For DNA quantitation, 260-, 280- and 340 nm measurements are all captured with a single flash.

## Advanced Alpha Technology detection

The combination of laser excitation and SDE detection gives the PHERAstar FSX unsurpassed flexibility and sensitivity for all AlphaScreen, AlphaLISA and AlphaPlex assays. The dedicated laser for Alpha Technology specifically excites donor beads at 680 nm and provides a broad dynamic range and an increased signal-to-noise ratio. In addition, the outstanding SDE detection system allows for reduced read times and higher sensitivities in all multiplex Alpha assays.




Assay-specific Optic Modules, fully-equipped, installed in seconds.



The cross-talk reduction package reduces unspecific light signals that are cause of higher background noise and variability.

## Cross-talk reduction


AlphaScreen/AlphaLISA and glow luminescence assays in 384- and 1536-well plates in particular are often negatively affected by stray light and cross-talk. Unwanted light may diffuse from adjacent wells to the detector, as well as through the plastic walls of the wells. If cross-talk is not reduced, low-signal wells might see more counts from nearby bright wells than from their actual signal. BMG LABTECH's cross-talk reduction package automatically applies an aperture to physically reduce non-specific signals diffusing from nearby wells, and mathematically corrects for light transmitted through the walls of a well.



Comparison of xenon lamp (10 flashes) and TRF laser (1 flash) for Signal-to-Blank (S/B), Z' value and read times with a TR-FRET assay.

## TRF-dedicated excitation laser

Dedicated lasers for excitation significantly improve assay performance and lower limits of detection as they yield a higher excitation energy at a specific wavelength. The PHERAstar FSX is equipped with a dedicated TRF laser for specific sample excitation at 337 nm. With 60 laser flashes per second, it allows for ultra-fast temporal resolution in TR-FRET measurements and even "flying mode" detection. For several applications, a single laser flash in flying mode provides enough energy to excite the donor molecules. Measurements can therefore take place without stopping plate movement, significantly reducing read times for an entire plate.



Unique TRF Decay Curve Monitoring feature for assay optimisation.

## Fine-tune your signal detection

Decay Curve Monitoring (DCM) is a unique feature of the PHERAstar FSX and a fundamental tool for assay development and performance fine-tuning in TRF, TR-FRET and Alpha Technology. Enabled by a dedicated photon counting detection system, DCM measures and displays the time-resolved emission curve of the fluorophore. The TR-FRET dedicated photon counting detection system enables simultaneous monitoring of both donor and acceptor decay curves with a time resolution of 2  $\mu$ s. Combined with the Integration Time Wizard, DCM helps to refine timing parameters, thus improving signal detection and reducing background noise.

## Advanced Well Scanning

Conventional readers excite and measure samples in the centre of the well. This can negatively affect data quality when measuring uneven cell layers or non-homogenous well samples. On the PHERAstar FSX, three different well scan modes enable robust data acquisition even from non-homogeneous samples such as adherent cells or bacteria. Orbital and spiral averaging automatically normalise for heterogeneous cell distributions. In these modes, the reader takes several measurements on a defined orbit and calculates an average for each well.


For higher resolution, matrix scan acquires up to 900 data points/well, displays each scan point graphically and creates a map for each well. Single scan points or entire sections can be easily removed upon detection.

## Unaffected by environmental changes


The new Advanced Assay Stability (AAS) system enables the PHERAstar FSX to match room temperature conditions and to be maintained at any temperature between 18 and 45°C. It achieves this by automatically heating or cooling the measurement chamber depending on the external environment. The AAS system provides a steady temperature that is unaffected by external environmental changes while the plate reader is operating, providing improved assay stability and more reliable data. The special design of the AAS system ensures that all plates in a screening campaign are measured under the same conditions. It provides a consistent environment for temperature-sensitive and kinetic assays, and reduces sample evaporation compared to conventional plate readers.

## Smart reagent injection

Reagent injection can be used to add a stimulus or inhibitor, and to initiate or stop kinetic and enzymatic reactions while detecting in real time the sample's reaction. On the PHERAstar FSX, two built-in reagent injectors allow for simultaneous detection and reagent injection. This ensures that even fast kinetic reactions are monitored without the loss of any data point. Injectors are readily accessible and housed within the instrument to protect any light-sensitive reagent. Extremely small dead volume and back flushing significantly reduce wastage of expensive reagent. The number of injections/well and injection volumes are adjustable for each well, allowing you to automatically produce dilution schemes and concentration gradients across the microplate.



Matrix scanning of GFP-expressing protoplasts. The scan shows how protoplasts (green colour) are not evenly distributed in the well.




The AAS system provides a stable temperature inside the PHERAstar FSX regardless of external temperature fluctuations (gray box).



The PHERAstar FSX can be equipped with a Stacker for mid-throughput applications.

## Data analysis made easy




MARS data analysis software for automated data reduction.

The software package includes the Reader Control and MARS data analysis interfaces. This multi-user software can be installed on as many computers as you require, without the need to purchase additional licenses.

The Reader Control software allows to define measurement protocols and acquire data. It is an extremely versatile interface for the straightforward execution of routine tasks, as well as the optimisation of complex operations.

MARS is designed to make data analysis simple and effective, and offers multiple data reduction possibilities such as:

- Standard Curve Wizard for a step-by-step standard curve calculation
- Automatic DNA/RNA concentration determination
- Data display as bar charts, box plots, violin plots etc.
- UV/vis spectral view and analysis
- Background and baseline correction
- Signal interpolation: linear or cubic spline
- Various curve fit models including linear, 4-/5-parameter, polynomial and user-defined fit models
- Enzyme kinetic analysis using various models
- EC<sub>50</sub> calculation with confidence intervals
- Binding rates and constants determination
- ANOVA, Student's t-test or multiple comparisons
- Performance evaluation: signal-to-blank, signal-to-noise, %CV, Z-prime, etc.
- Automatic data processing using predefined templates




NanoBRET competitive binding experiments of propranolol-BY630 with increasing concentrations of known unlabelled  $\beta$ 2AR ligands.

The software package comes with flexible data export (Excel, ASCII) and integration capabilities, and is fully compliant with FDA regulation 21 CFR Part 11.

## Applications hub

A perfectly engineered instrument is only part of the solution, it needs to effectively perform all of the leading applications. We continuously work with all major reagent companies to develop protocols and improve instrument settings for their existing assays and their newest kits. The PHERAstar FSX is a user-friendly and flexible instrument that provides excellent performance in all applications, including NanoBRET™, AlphaLISA® and TR-FRET assays, as shown in the examples. Applications include:

- Protein-protein interactions
- Kinetic binding assays
- Compound and inhibitor screening



Proof of performance of an AlphaLISA-based SARS-CoV-2 Spike Protein - ACE2 competitive binding assay.

- DNA, RNA, and protein quantification
- Enzyme activity and kinetic assays
- Cell-based assays
- And much more ...

Our comprehensive online application database reflects more than 30 years of expertise and innovations. Over 8,000 published entries of peer-reviewed articles and application notes demonstrate the flexibility and versatility of our readers, and their use in chemical and biological sciences.


## Automation

The PHERAstar FSX offers excellent robotic integration capabilities, multi-user control, digital signature and FDA 21 CFR part 11 compliance. Its robotic software interface make it easy to integrate into all leading robotic platforms. The instrument comes with three integrated microplate barcode readers and can be equipped with a Stacker.

## Support and training

BMG LABTECH operates globally through an extensive network of subsidiaries and trained distributors. Customers can rely on qualified support and assistance with regard to software, assay development, or general enquiries related to the PHERAstar FSX and all our other microplate reading solutions.

AlphaTechnology includes AlphaScreen, AlphaLISA, and AlphaPlex. These are registered trademarks of PerkinElmer, Inc.  
 HTRF is a registered trademark of Cisbio Bioassays.  
 LanthaScreen is a registered trademark of Thermo Fisher Scientific.  
 Transcreener is a registered trademark of Bellbrook Labs.  
 DLR and NanoBRET are trademarks of Promega Corp.  
 Mycoalert is a registered trademark of Lonza.  
 THUNDER is a trademark of Bioauxilium.



TR-FRET comparison between the PHERAstar FSX and another HTS microplate reader, using a THUNDER human IL-1β standard curve.



The PHERAstar FSX can include all or any combination of features/options/accessories listed below at purchase. Upgrading with additional features/options/accessories may be possible after purchase. Contact your local representative for more details or a quote.



**Headquarters Germany**  
**BMG LABTECH GmbH**  
 Allmendgrün 8  
 77799 Ortenberg  
 Tel. +49 781 96968 -0  
 sales@bmglabtech.com

**Australia**  
 BMG LABTECH Pty. Ltd.  
 2/24 Carbine Way  
 Mornington, Victoria, 3931  
 Tel. +61 3 5973 4744  
 australia@bmglabtech.com


**France**  
 BMG LABTECH SARL  
 7, Rue Roland Martin  
 94500 Champigny s/Marne  
 Tel. +33 1 48 86 20 20  
 france@bmglabtech.com

**Japan**  
 BMG LABTECH JAPAN Ltd.  
 1-6-2, Shimo-cho  
 Omiya-ku  
 330-0844 Saitama City  
 Tel. +81 48 647 7217  
 japan@bmglabtech.com

**UK**  
 BMG LABTECH Ltd.  
 8 Bell Business Park  
 Smeaton Close  
 Aylesbury, Bucks  
 HP19 8JR  
 Tel. +44 1296 336650  
 uksales@bmglabtech.com

**USA**  
 BMG LABTECH Inc.  
 13000 Weston Parkway  
 Suite 109  
 Cary, NC 27513  
 Tel. +1 877 264 5227  
 usa@bmglabtech.com

[www.bmglabtech.com](http://www.bmglabtech.com)



Detection modes	UV/vis absorbance Fluorescence intensity - incl. FRET Luminescence (flash and glow) - incl. BRET Fluorescence polarization Time-resolved fluorescence TR-FRET AlphaScreen®/AlphaLISA®/AlphaPlex™	
Measurement modes	Top and bottom reading Endpoint and kinetic Sequential multi-excitation Sequential multi-emission Simultaneous Dual Emission Spectral scanning (absorbance) Real-time ratiometric measurements Well scanning	
Microplate formats	Up to 3456-well plates, user-definable LVis Plate with 16 low-volume microspots (2 µL)	
Microplate carrier	Robot compatible	
Light sources	High energy xenon flash lamp Dedicated laser for AlphaScreen®/AlphaLISA®/AlphaPlex™ Dedicated laser for TRF and TR-FRET	
Detectors	Two matched pairs of photomultiplier tubes (PMTs), dedicated to different detection modes CCD spectrometer	
Wavelength selection	Optic Modules: up to six application-specific modules (includes all filters, dichroics, etc., for an application) For top and bottom reading UV/vis absorbance spectrometer: full spectrum or up to 8 distinct wavelengths in < 1 sec / well	
Optical path	Top and bottom: free-air optical light path guided by motor-driven mirrors and dichroics	
Z-adjustment	Automatic focal height adjustment (0.1 mm resolution)	
Spectral range	<b>Filters</b>	230 - 750 nm or 230 - 900 nm for FI, FP 230 - 750 nm for LUM 230 - 900 nm for TRF
	<b>Spectrometer</b>	220 - 1000 nm for ABS; wavelength precision: < 0.5 nm
Sensitivity	<b>FI (top)</b>	< 0.15 pM FITC (384sv black, 20 µL) < 0.5 pM FITC (1536 black, 8 µL)
	<b>FI (bottom)</b>	< 1.0 pM FITC (384g black, 50 µL)
	<b>FP</b>	< 0.5 mP SD at 1 nM FITC (384sv black, 20 µL) < 1.5 mP SD at 1 nM FITC (1536 black, 8 µL)
	<b>TRF with laser</b>	< 5 fM europium (384, 80 µL) < 15 fM europium (1536, 8 µL)
	<b>HTRF® (black and white microplates)</b>	Reader Control Kit [Eu] after 18h (384sv, 20 µL) Delta F > 1100 % (High Calibrator) Delta F > 25 % (Low Calibrator)
	<b>LUM</b>	< 0.4 pM [ $< 8$ amol/well ATP, 384sv white, 20 µL] Dynamic range: 9 decades
	<b>AlphaScreen® with laser</b>	< 5 pM [ $< 100$ amol/well P-Tyr-100, 384sv white, 20 µL]*
	<b>ABS with spectrometer</b>	Selectable spectral resolution: 1, 2, 5, and 10 nm OD range: 0 to 4 OD; photometric resolution: 0.001 OD Accuracy: < 1% at 2 OD Precision: < 0.5% at 1 OD and < 0.8% at 2 OD Linearity: < 0.8% at 2.0 OD
Read times	<b>Flying mode (1 flash)</b>	14 sec (384) 27 sec (1536)
	<b>10 flashes</b>	38 sec (384) 1 min 52 sec (1536)
	<b>50 flashes</b>	1 min 29 sec (384) 5 min 18 sec (1536)
Reagent injection	Up to 2 built-in reagent injectors Injection at measurement position (6 to 384-well) Individual injection volumes for each well (3 to 500 µL, optionally up to 2 mL) Variable injection speed up to 420 µL/sec Up to four injection events per well Reagent back flushing	
Shaking	Linear, orbital, and double-orbital with user-definable time and speed	
Integrated barcode reader	Up to three integrated barcode readers	
Incubation	+5°C above ambient up to 45°C. The upper heating plate of the incubation chamber operates at +0.5°C than the lower plate. This prevents condensation build-up on the lid or sealer With Advanced Assay Stability (AAS) system: from 18 - 45°C	
Software	Multi-user Reader Control and MARS data analysis software included FDA 21 CFR Part 11 compliant	
Dimensions	Width: 45 cm, depth: 51 cm, height: 47 cm; weight: 49 kg With AAS system: width: 48 cm, depth: 68 cm, height: 47 cm; weight: 62 kg	
	<b>Accessories</b>	
Stacker	Plate handler for up to 50 microplates - continuous loading feature	
ThermoStar	Microplate incubator and shaker	
LVis Plate	Microplate designed to measure 16 low-volume (2 µL) samples and standard cuvettes. Incorporates NIST-traceable filters and holmium oxide standards for instrument performance test. Sensitivity: 2 ng/µL dsDNA	
Upgrades	Please contact your local representative for upgrades including options such as detection modes, reagent injectors, etc.	

\* Limit of detection: < 100 amol of biotinylated and phosphorylated polypeptide (P-Tyr-100 assay kit, PerkinElmer, #6760620C), measured in white 384 small volume microplates

Limit of detection (sensitivity) was calculated according to the IUPAC standard:  $3x(SD_{blank})/slope$   
 Specifications are subject to change without notice.

All mentioned products are for research use only.

© 2023 All rights reserved. All logos and trademarks are the property of BMG LABTECH.