

Quantitative Real-time PCR system

(ZIP-96V)



Molecular

www.zybio.com

Quantitative Real-time PCR system (ZIP-96V)

-

ZIP-96V

Quantitative real-time PCR system ZIP-96V is 96-well real-time PCR instrument based on Fresnel lens optical signal acquisition technology, time-resolved signal separation technology and unique temperature control technology. It is designed for multiple fields of basic research in medicine and biology, such as pathogen detection, genotyping, gene therapy drug, gene expression, food safety testing, public health, animal health, etc.

Features



• Six independent temperature control blocks improve the accuracy and uniformity of temperature control

 Unique time-resolved signal separation technology, no multi-color crosstalk

signal sensitivity and signal to noise ratio

Multiple Applications

- · Basic scientific research
- Pathogen detection
- Gene expression
- Food safety testing
- GenotypingPublic health
- Gene therapy drug
- Animal health

Four fluorescence Channels, support most of common dyes

Detection Reporter dye:



Technical parameters

Basic parameters

Model	ZIP-96V		
Dimensions(W×D×H, mm)	425*320*205		
Weight	13.5Kg		
Throughput	Up to 96 tests per run		
Reaction volume	20-100µL		
	0.2ml PCR tube/		
Sample Format	8-strip PCR tube/ 96-well plate		
Supported dye	Channel 1: FAM		
	Channel 2: VIC, HEX		
	Channel 3: ROX		
	Channel 4: CY5		
Fluorescent channels	4 detection channels		
Thermal cycle technology			
Temperature control	Peltier		
Temperature control block	6 blocks		
Temperature range	4~99°C		
Average heating rate	> 2.5°C/s		
Average cooling rate	> 2.0°C/s		
Temperature fluctuation	±1°C		
Temperature accuracy	≤ ±0.5°C		
Hot-lid temperature range	30~110°C		

Optical system

High efficiency LEDs		
Photodiode		
460nm~480nm, 525nm~545nm,		
575nm~595nm , 600nm~620nm		
505nm~525nm, 555nm~575nm,		
600nm~620nm, 660nm~680nm		
< 2%		
r > 0.99		
Qualitative/		
Absolute Quantification Analysis		
All-channel scanning		
Excel		
Multiple printing templates/		
customized template		
100-240 V AC, 50-60 Hz		
Windows 7/10		
USB		
10-30°C		
20-85%		

Classic Examples

SARS-CoV-2 Nucleic Acid Detection (PCR-Fluorescent Probe Method) qPCR technology was used to detect and diagnose suspected case of COVID-19.

Detection kit:

SARS-CoV-2 nucleic acid detection kit (Zybio)

· Isolation kit and system:

nucleic acid extraction kit (Zybio), nucleic acid isolation system EXM3000 (Zybio)

Assay method:

Take 200µL of sample (throat swab preservation solution) to be tested, and extract RNA from the sample with nucleic acid extraction kit and isolation system EXM3000. According to the instruction of nucleic acid detection kit, mix extracted RNA, the positive control product, and the negative control product with PCR reaction solution and enzyme respectively, mix well and centrifuge and then do amplification and detection.

. Amplification parameters:

Steps		Temperature	Time	Cycle
1	UNG reaction	37°C	1 min	1
2	Reverse transcription	50°C	5 min	1
3	Initial denaturation	95°C	2 min	1
4	Denaturation	95°C	5 sec	45
5	Amplification and fluorescence detection	O°C	30 sec	

Fluorescence Detection: Step 5.

Report Fluorescence Setting: FAM, ROX, VIC;

Quenching Fluorescence Setting: None; Passive Reference Setting: None.

· Cut-off Value:

result is clamed positive when the Ct of two targets (FAM, ROX) <40.

· Result analysis:

Positive control and negative control are normal, sample to be tested is positive. The corresponding amplification curves are as follows:





Zybio Inc.

Address: Floor 1 to Floor 4, Building 30, No.6 of Taikang Road, Block C of Jianqiao Industrial Park, Dadukou District, Chongqing, China 400082 Tel: +86-23 6865 5509 Fax: +86-23 6869 9779 Email: info@zybio.com Website: www.zybio.com

EN-C-FZ-ZIP-96V-20210105H

Molecular

www.zybio.com